## HEINE EL3 LED Examination Light

**HEINE Optotechnik GmbH & Co. KG**  
Kientalstr. 7 · 82211 Herrsching · Germany  
Tel. +49(0)8152/38-0  
Fax +49(0)8152/38-202  
E-Mail: [info@heine.com](mailto:info@heine.com) · [www.heine.com](http://www.heine.com)  
med 0315    2016-05-23

# HEINE EL3 LED Examination Light

DEUTSCH	3
ENGLISH	5
FRANÇAIS	7
ESPAÑOL	9
ITALIANO	11
SVENSKA	13
NEDERLANDS	15
DANSK	17
NORSK	19
SUOMI	21
PORTUGUÊS	23

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# HEINE EL3 LED Examination Light

 Please read and follow these instructions for use of and keep them for future reference.

## Intended Use

The EL3 LED is a mains operated medical examination light for illuminating the body's surfaces and cavities during medical examinations in closed examination rooms.

 The device is intended only for medical examinations in which failure of the light will not pose a risk to the patient. A main source of lighting must always be switched on.

## For U.S. only:

 Federal law restricts this device to sale by or on the order of a Physician or Practitioner.

## Warnings and Safety Information

 **CAUTION!** Indicates potential hazardous situations. Ignoring the corresponding instructions may lead to dangerous situations of mild to moderate extent. (Background color yellow; foreground color black).

 **NOTE!** Note indicates valuable advice in terms of installation, operation, maintenance or repair. Notes are important, but not related to hazardous situations.

## Product overview

Gooseneck:



- 1 Gooseneck (flexible component)
- 2 Gooseneck (rigid component)
- 3 Gooseneck (flexible component)

Lamp head:

- 4 Hand grip


Lamp base:

- 5 ON/OFF switch (system)
- 6 Connection socket for the power cable
- 7 Installation cover
- 8 Connection for the mounting systems

For mounting the mains adaptor to the wall, clamp or stand mount, push up the cover (7) and take it off. Use the screws supplied with the product to fasten the examination light to the respective mounting system (clamp or stand mount) or use the wall plugs supplied to mount the lamp to the wall.

The screws for the mounting systems are located in the corners of the mounting unit.

### ① Clamp mount:

Please use the screws enclosed with the product to attach the clamp mount to the base of the lamp. This can be fastened to an even surface (e.g. table top) with a maximum thickness of 40 mm or to a circular tube with a maximum diameter of 30 mm. Please check that the mounted lamp is securely in place before using the lamp.

### ② Stand mount:

This mount system is designed for tubes with a diameter of 25 mm. Please use the screws enclosed with the product to attach the stand mount to the base of the lamp.

### ③ Wall mounting:

The rubber mat can be used as a drilling template for marking the holes on the wall and is intended to level out uneven areas on the wall for the mains adaptor.

Please use the screws in combination with the wall plugs supplied with the product when mounting the device to a suitable wall.

## Setting up

 The EL3 LED must not be used for eye examinations or as a surgical light.

The EL3 LED is intended for use in closed spaces. It is not designed for mobile use or for use outdoors. Improper use can result in damage to the device or pose a risk to the user.

Prior to each use, check that the EL3 LED is functioning properly as intended. Do not use the device if visible damage is apparent.

The device does not have an application part and therefore care must be taken to ensure that the examination light does not come into contact with the patient!

The EL3 LED must not be used in strong magnetic fields such as within the vicinity of MRI equipment.

The EL3 LED must not be used in oxygen-enriched environments.

The EL3 LED must not be used in an environment of flammable anaesthetic agents (class AP) or with oxidising agents (class APG).

Keep a minimum distance of 20 cm between the examination light and the examined object (patient).

Do not cover the examination light.

Ensure that the device is handled properly and carefully.

Do not let the power cable become taut as this could damage the device or pose a tripping hazard.

Only after the device has been securely mounted, should the power cable supplied with the device be connected by plugging it into the mains connection socket (6). To use the device, switch it on by pressing the power switch (5) so that the switch is in the ON position I and the switch lights up green.

## Installation instructions

When mounting the accessories, please follow the enclosed mounting instructions for the mounting systems.

 To avoid risk of electric shock, this equipment must only be connected to a supply mains with protective earth.

To ensure all-pole disconnection from the mains at all times, the ME device must be installed so that the power cord is accessible and disconnectable.

When using the examination light with a stand, clamp or wall mount, ensure the examination light is standing or stably attached. Check that the attached examination light is stable and secure before each use.

Use the accessories supplied with the examination light to ensure the examination light is held securely and properly in the respective mount.

When moving the lamp, ensure that the power cable is not taut. Unplug the power cable when moving the lamp over longer distances.


When mounting the lamp to the wall, ensure that the subsurface is able to bear the mechanical stress.

The wall plugs supplied with the product are universal wall plugs and suitable for most building materials (e.g. concrete, solid brick, plasterboard or gypsum cardboard). Before mounting the device, check whether you need to use a special wall plug for the subsurface of your wall.

Before mounting the device to the wall, check that there are no cables or pipes running under the planned drill holes that could get damaged.

## Positioning

- ⚠ To avoid contamination, position the examination light prior to the examination or by an assistant during the examination.

The swivel arm of the EL3 LED examination light has two flexible (1, 3) and one rigid section (2).

To position the light precisely, hold the light head by the hand grip (4) and the gooseneck by the flexible part (1) and bend the gooseneck into the desired position. Fine adjustments can be made via the hand grip (4).

Do not bend the gooseneck into a <math><90^\circ</math> degree angle.

The picture shows the positioning of the lamp during transport on a stand.

## Operation

- ⚠ Keep a minimum distance of 20 cm between the examination light and the examined object (patient). The optimal working distance is 30cm.

## Hygienic Reprocessing

Instructions on hygienic reprocessing must be adhered to, based on national standards, laws and guidelines.

**Classification according to KRINKO:** non-critical

**Spaulding Classification USA:** non-critical

- ⚠ Allow the device to cool down before reprocessing.

Before cleaning, disconnect the device from the power source.

In the event of suspected contamination, carry out hygienic preparation of the instrument.

The described cleaning and disinfection measures do not replace the specific rules applicable for the establishment.

HEINE Optotechnik only approves the resources and procedures named in these instructions for use.

Cleaning and disinfection may only be carried out by personnel with sufficient hygienic knowledge.

Observe the instructions of the manufacturer of the reprocessing media.

Do not use any spray disinfectant, or any dripping-wet or heavily-foaming cloths. Do not prepare mechanically.

## Procedure

Clean and disinfect the EL3 LED manually (wipe clean and wipe disinfect).

Recommended agents:

Alcoholic (eg Incides<sup>®</sup> tissues), quaternary ammonium compounds (eg Microbac<sup>®</sup> tissues) or ortho-phthalaldehyde (OPA example Cidex<sup>®</sup>)

Conditioning residues on optics may be removed using a low-bacterial, dry, micro-fibre cloth.

## Changing the light source

- 🔧 The LED unit cannot be replaced. In case of a fault or defect, please contact HEINE.


## Replacing fuses

- ⚠ Switch off the examination light and remove the power cable from the device. The fuses are located in the fuse holder of the panel mount connector. Use only the fuses specified in the technical data as replacement fuses.

## Technical safety inspection

- ⚠ The power supply unit and power cable should be inspected for damage and that they are in good working order at regular intervals.

The technical safety inspection must be carried out every two years in accordance with the generally accepted rules of technology.

- ① For the electrical inspection, the installation cover must be removed first.
- ② The testing point is located on the left hand side and is labelled accordingly.

## Maintenance and Service

The instrument does not require maintenance or service.

The device may not be modified in any way.

For repair, send the device to HEINE or an authorized service partner only, to avoid the risk of electrical shocks, burns or injuries due to improper repair work or handling.

## General Warnings

- ⚠ Check the correct operation of the device before use! Do not use the device if there are visible signs of damage.

Do not modify the device.

Use only original HEINE parts, spare parts, accessories and power sources.

Repairs shall only be carried out by qualified persons.

## General Notes

- 🔧 The warranty for the entire product is invalidated if non-genuine HEINE products or non-original parts are used and if repairs or modifications are made to the device by persons not authorized by HEINE. For more information, please visit [www.heine.com](http://www.heine.com).

Condensation can form when moving from a cold to a warm environment. Wait until the unit has warmed up and is absolutely dry before using it.

## Disposal

- 🗑 The product must be recycled as separated electrical and electronic devices. Please observe the relevant state-specific disposal regulations.

## Electromagnetic Compatibility

Medical electric devices are subject to special precautionary measures with regard to electromagnetic compatibility (EMC). Portable and mobile high frequency communication equipment can affect medical electric devices.

- ⚠ This is a device in the domestic environment, this device may cause radio interference, so that it may be necessary in this case, to take appropriate remedial measures, as e.g. orientation, new arrangement or shielding of the device or restrict the connection to the site.

The use of accessories, converters or cables other than the ones specified by HEINE might lead to increased emission reduced electrical immunity of the medical equipment.

The device may not be stacked directly near or used directly beside other devices. If the device is to be operated in a stack or with other devices, the device should be watched to ensure it operates properly in this location.

The appendix contains following tables:

- Guidance and Manufacturer's declaration – Electromagnetic immunity
- Technical specification
- Explanation of the used symbols

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**Guidance and manufacturer's declaration – electromagnetic emissions**

The device is intended for use in the electromagnetic environment specified below.  
The customer or the user of the device should assure that it is used in such environment.

Emission test	Compliance	Electromagnetic environment – Guidelines
RF emissions CISPR11	Group 1	The device uses RF energy only for its internal function. Therefore, RF-emission is very low and it is unlikely that any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic Emissions IEC 61000-3-2	Class A	Symmetrical three-phase devices and other devices.
Voltage Fluctuations/ Flicker Emissions IEC 61000-3-3	Passed	

**Guidance and manufacturer declaration - Electromagnetic immunity**

The device is intended for use in the electromagnetic environment specified below.  
The customer or the user of the device should assure that it is used in such environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – Guideline
Electrostatic discharge (ESD) acc. to IEC 61000-4-2	± 6 kV contact discharge ± 8 kV air discharge	± 6 kV contact discharge ± 8 kV air discharge	Floors should be wood, concrete or covered with ceramic tiles. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for mains cables ± 1 kV for input and output lines	± 2 kV for mains cables ± 1 kV for input and output lines	The supply voltage quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV voltage phase – phase, ± 2 kV voltage phase – earth	± 1 kV voltage phase – phase ± 2 kV voltage phase – earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	< 5% UT, (>95% dip in UT) for 1/2 period 40% UT, (60% dip in UT) for 5 periods 70% UT, (30% dip in UT) for 25 periods <5% UT, (>95% dip in UT) for 5 seconds	< 5% UT, (>95% dip in UT) for 1/2 period 40% UT, (60% dip in UT) for 5 periods 70% UT, (30% dip in UT) for 25 periods <5% UT, (>95% dip in UT) for 5 seconds	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered by a UPS (uninterruptible power supply) or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	10 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Comment: UT is the a.c. supply voltage prior to application of the test level.			

**Guidance and manufacturer's declaration – electromagnetic immunity**

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment Guidelines
Conducted RF IEC 61000-4-6	3 Veff 150 kHz to 80 MHz	10 V eff	Portable and mobile RF communication equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Radiated HF IEC 61000-4-3	3 V/m 80MHz to 2,5GHz	10 V/m	<p><b>Recommended separation distance:</b></p> $d = 3,5/3 \cdot \sqrt{P/W}$ $d = 3,5/3 \cdot \sqrt{P/W} \text{ 80 MHz to 800 MHz}$ $d = 7/3 \cdot \sqrt{P/W} \text{ 800 MHz to 2,5 GHz}$ <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey<sup>a</sup>, should be less than the compliance level in each frequency range.<sup>b</sup> Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

Note 1: At 80Hz and 800MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

**a** Field strengths from fixed transmitters, such as base stations for radio (cellular / cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy.  
To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation.  
If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.

**b** Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V / m.

### Recommended separation distances for portable and mobile RF communication equipment and the device

The device is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communication equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d = 3,5/3 \cdot \text{SQRT}(P)$	80 MHz to 800 MHz $d = 3,5/3 \cdot \text{SQRT}(P)$	800 MHz to 2,5 GHz $d = 7/3 \cdot \text{SQRT}(P)$
0.01	0.1	0.1	0.2
0.1	0.4	0.4	0.7
1	1.2	1.2	2.3
10	3.7	3.7	7.4
100	11.7	11.7	23.3

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in metres (m) can be determined using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

### Technical specification

Mains power	100–240 VAC, 50/60 Hz
Power consumption	25 VA
Fuse sets	2 x T 1.25AH 250V
Average LED operating life	typ. 30,000 hours
Minimum distance from light source to patient	20 cm
Working distance from light source to patient	30 cm
Central illuminance $E_c$	typ. 52,700lx*
Total irradiance $E_e$	typ. 216 W/m <sup>2</sup> *
Colour temperature	5000 K $\pm$ 500 K
Colour rendering index Ra	$\geq$ 90
Special colour index R9	typ. 90
Environmental conditions for operation	+10°C to +35°C 30% to 75% rel. humidity 700hPa to 1060hPa
Environmental conditions for storage	+5°C to +45°C 45% to 80% rel. humidity 500hPa to 1060hPa
Environmental conditions for transport	-20°C to +50°C 45% to 80% rel. humidity 500hPa to 1060hPa
Protection class	I
Operation altitude	$\leq$ 3000 m
Protection provided by enclosures (IP code)	IP20
Risk group (according to DIN EN 62471)	1 (low risk)
Weight incl. wall mount	approx. 1.8 kg
Weight of clamp mount	approx. 0.2 kg
Length of flexible arm	1 m

\*at a working distance of 30 cm, depending on  $\theta$

UL60601-1 / CAN/CSA C22.2No.601.1  
ANSI/AAMI ES60601-1 (2005)  
CAN/CSA-C22.2 No. 60601-1 (2008)  
28KY

MEDICAL — GENERAL MEDICAL EQUIPMENT  
AS TO ELECTRICAL SHOCK, FIRE AND MECHANICAL HAZARDS ONLY  
IN ACCORDANCE WITH UL60601-1 / CAN/CSA C22.2No.601.1  
ANSI/AAMI ES60601-1 (2005) CAN/CSA-C22.2 No. 60601-1 (2008)  
28KY

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#### Explanation of utilized symbols

The following symbols are used on the device or on the packaging:

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