HEINE QUALITY MADE IN GERMANY	DEUTSCH
	ENGLISH
	FRANÇAIS
	ESPAÑOL
	ITALIANO

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HEINE Ootoscope G100/G100 LED Slit Illumination Head

! Please read and follow these instructions for use of and keep them for future reference.

Intended Use

The HEINE Ootoscope G100/G100 LED Slit Illumination Head is intended for the examination of the ear, the external auditory canal and the tympanum. The instrument has a swivel lens for optical magnification and an internal powered illumination unit. The device is intended for short term treatment and can also be used for instrumentation with non-active medical devices. The instrument may only be used by qualified medical personnel and in professional healthcare facilities.

! **For U.S. only:** **Federal law restricts this device to sale by or on the order of a Physician or Practitioner.**

Warnings and Safety Information

! **CAUTION!** Indicates potential hazardous situations. Ignoring the corresponding instructions may lead to dangerous situations of mild to moderate extent. (Background color yellow; foreground color black).

! **NOTE!** Note indicates valuable advice in terms of installation, operation, maintenance or repair. Notes are important, but not related to hazardous situations.

Product overview (see figure)

1 Slit Illumination Head	6 Spreadable speculum
2 HEINE Handle	7 Short, closed specula
3 Swivel lens	8 HEINE UniSpec® Disposable Specula
4 Long, closed metal specula	9 Adaptor for UniSpec Specula
5 Slit metal specula, open	

Setting up

To set up the HEINE device, plug it onto the HEINE power source.

Verify that the lamp voltage complies with the supply voltage of the power source.

The coloured marking on the bottom of the lamp shows you the lamp voltage:

White ring = HEINE XHL® 2.5 V bulb

only for use with the HEINE BETA Battery handle

Red ring = HEINE XHL® 3.5 V bulb

Black ring = HEINE LED illumination 3.5 V

Only for use with the HEINE BETA rechargeable handle and other HEINE power sources.

Operation

The Slit Illumination Head may only be inserted into the auditory canal when a tip or specula is fitted.

Place the HEINE reusable specula or specula in such a way onto the Slit Illumination Head

that the guide pin fits into the bayonet slit. Tighten the tip by twisting it slightly to the right.

HEINE UniSpec® disposable specula can be attached by using an adaptor.

Swivel Lens

The instrument head is fitted with a swivel lens of approx. 1.6 x magnification. The Slit Illumination Head may be used with all specula, except the 150 mm long, closed metal specula, as this speculum may affect the focusing.

Application duration

The Slit Illumination Head is intended for short-term treatment.

XHL: Application duration < 2 minutes with a break between applications of 15 minutes.

LED: Application duration < 6 minutes with a break between applications of 10 minues.

The setting up and operation of the HEINE power supplies are described in a separate instruction of use.

Hygienic Reprocessing

Instructions on hygienic reprocessing must be adhered to, based on national standards, laws and guidelines.

Classification according to KRINKO: noncritical

Spaulding Classification USA: noncritical

! Lens and illumination unit must be disassembled from the instrument head before hygienic reprocessing. These cannot be automatically reprocessed or sterilized.

Reusable specula and specula are to be reprocessed hygienically after every application.

In the event of suspected contamination of the Slit Illumination Head, carry out hygienic reprocessing of the instrument.

Sterilisation and automated reprocessing of the Slit Illumination Head should only be carried out after the treatment of high-risk patients or for semicritical applications.

The described cleaning and disinfection measures do not replace the specific rules applicable for the establishment.

HEINE Optotechnik only approves the reprocessing agents and procedures named in these instructions for use.

Cleaning and disinfection may only be carried out by personnel with sufficient hygienic knowledge.

Observe the instructions of the manufacturer of the reprocessing media.

Do not use spray disinfection and ultrasonic cleaning.

The sterilization of the slit illumination head via moist heat and mechanical reprocessing above 65 °C can lead to a reduction of the life expectancy of the fiber optic guide and possible losses in light transmission.

HEINE UniSpec® Disposable Specula are for single use only.

Procedure

Slit Illumination Head for noncritical applications

Cleaning and disinfection

Clean and disinfect the Slit Illumination Head and the loupe with a damp cloth (wipe cleaning and wipe disinfection).

Recommended agents

Cleaning agent: enzymatic (e.g. neodisher® MediClean)

Disinfectant agent: alcoholic (e.g. Incides® N), or quaternary ammonium compounds (e.g. Mikrobac® tissues)

Slit Illumination Head for semicritical applications as well as reusable specula and specula

Cleaning and disinfection

Clean and disinfect the Slit Illumination Head, the reusable specula and specula manually (in an immersion bath) or automic (Method according to EN ISO 15883).

Recommended agents

Cleaning agent: enzymatic (e.g. neodisher® MediClean)
Disinfectant agent: alcoholic (e.g.rotasept®), agent ortho-phthalaldehyde (e.g. Cidex® OPA) or quaternary ammonium compounds (e.g. neodisher® Septo MED)

Sterilization

Recommended sterilization programs

Humid heat (ISO 17665)

132 – 134 °C, 3 min

Fractionated vacuum procedure or gravity displacement procedure

Changing the light source

! Allow the device to cool down before changing the bulb.

XHL and LED bulb may not be interexchanged, as the corresponding instrument heads have different dimensions.

- Remove the instrument head from the handle.

XHL: Remove the XHL lamp by pulling on the narrow ridge of the lamp. Insert the new lamp into the guide tube as far as it will go.

LED: Use the supplied screwdriver to loosen the small screw in the connector and carefully pull out the LED lamp by the contact pin. Insert the new LED as far as it will go and secure it by lightly screwing in the small screw. Do not overtighten the screw.

Maintenance and Service

The instruments do not require maintenance or service.

We as manufacturer do not stipulate a safety inspection in accordance with § 6 Medical Device Operator Ordinance (MPBetreibV), ref. attachment 1, for the HEINE product G100/G100 LED. The expected operational life of the HEINE Ootoscopes is 7 years.

General Warnings

! Check the correct operation of the device before use! Do not use the device if there are visible signs of damage or the light begins to flash.

Do not use the device in fire- or explosive risk area (e.g. oxygen saturated or anaesthetic environments).

This product is not allowed to enter or be used in areas with strong magnetic fields e.g. MRI scanners.

Do not modify the device.

Use only original HEINE parts, spare parts, accessories and power sources.

Repairs shall only be carried out by qualified persons.

Do not shine light of the instrument directly into the eyes.

Choose a tip/specula size suitable for the patient.

Check the reusable specula and speculum before each use for rough surfaces and edges.

In order not to exceed the surface temperature of 41 °C of the application part, keep within the application duration.

Do not separate the diagnosis head from the handle in the patient environment.

General Notes

! The warranty for the entire product is invalidated if non-genuine HEINE products or non-original parts are used and if repairs or modifications are made to the device by persons not authorized by HEINE. For more information, please visit www.heine.com.

Disposal

! The product must be recycled as separated electrical and electronic devices. Please observe the relevant state-specific disposal regulations.

Electromagnetic compatibility

Medical electric devices are subject to special precautionary measures with regard to electromagnetic compatibility (EMC). Portable and mobile high frequency communication equipment can affect medical electric devices.

! This is a device in the domestic environment, this device may cause radio interference, so that it may be necessary in this case, to take appropriate remedial measures, as e.g. orientation, new arrangement or shielding of the device or restrict the connection to the site. The use of accessories, converters or cables other than the ones specified by HEINE might lead to increased emission and reduced electrical immunity of the medical equipment.

The device may not be stacked directly near or used directly beside other devices. If the device is to be operated in a stack or with other devices, the device should be watched to ensure it operates properly in this location.

The appendix contains following tables

- Guidance and manufacturer's declaration – Electromagnetic immunity

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